

# User Manual

Release 1.11

## PAT-Cell-Opto-10

Electrochemical test cell



The information in this manual has been carefully checked and believed to be accurate; however, no responsibility is assumed for inaccuracies.

EL-Cell GmbH maintains the right to make changes without further notice to products described in this manual to improve reliability, function, or design. EL-Cell GmbH does not assume any liability arising from the use or application of this product.

**EL-Cell GmbH**

Tempowerkring 8

21079 Hamburg - Germany

**phone:** +49 40 79012-737

**fax:** +49 40 79012-736


**e-mail:** [support@el-cell.com](mailto:support@el-cell.com)

**web:** [www.el-cell.com](http://www.el-cell.com)

# Content

<b>1 Product Description .....</b>	<b>4</b>
<b>2 Technical data.....</b>	<b>6</b>
<b>3 Intended use.....</b>	<b>6</b>
<b>4 Safety Precautions .....</b>	<b>6</b>
<b>5 Assembly .....</b>	<b>7</b>
<b>6 Connecting the test cell .....</b>	<b>16</b>
<b>7 Disassembly and Cleaning.....</b>	<b>17</b>
<b>8 Unpacking.....</b>	<b>17</b>
<b>9 Components .....</b>	<b>19</b>
<b>10 Technical Support .....</b>	<b>21</b>
<b>11 Warranty.....</b>	<b>21</b>

# 1 Product Description




The PATCell-Opto-10 test cell is designed for operando characterization of electrodes using optical methods such as light microscopy or Raman spectroscopy in reflection mode.

The test cell utilizes the cableless PAT socket for cell connection. This way, it can be directly plugged into a PAT battery tester like the PAT-Tester-x-8 or a PAT docking station to connect it to a third-party battery tester.

## Features


- 2- and 3-electrode cell with optical window for aprotic and aqueous electrochemistry.
- High cycling stability due to glass-to-metal seals
- Suited for light microscopy and Raman spectroscopy
- Cableless cell connection via PAT socket

Working principle of the PAT-Cell-Opto-10 (side-by-side setup of electrodes):



## 2 Technical data

- Length: 55 mm
- Width: 55 mm
- Height: 32.5 mm
- Weight: 0.3 kg
- Operational temperature:  $-20^{\circ}$  to  $+70^{\circ}\text{C}$
- Electrode R diameter:  $\leq 10$  mm
- Dead volume:  $1.8\text{ cm}^3$



## 3 Intended use

The PAT-Cell-Opto-10 test cell is an electrochemical measuring instrument designed for use in a laboratory environment. It may only be used by trained personnel and only as described in this manual.

## 4 Safety Precautions

Use proper safety precautions when using hazardous electrode materials and electrolytes. Wear protective glasses and gloves to protect you against electrolytes that may accidentally spill out during filling and disassembly. Upon cell disassembly, dispose of all materials properly. Metallic lithium and some insertion compounds may decompose heavily in contact with water and other solvents and can cause a fire.

## 5 Assembly

**Please note:** All assembly steps are to be carried out in an inert glove box atmosphere. All components used are to be dried upfront in a vacuum oven at 80°C (120°C for parts made of PPS) for at least 12 hours.

Once fully assembled, the cell is hermetically sealed so that it may be operated in an ambient atmosphere.

The test cell can be used in several different configurations, which mainly differ in the type and size of the working electrode used, the charging geometry (the position of the working and counter electrode relative to each other), and the connection to the potentiostat (2- or 3-electrode connection). In the following, the cell assembly is described in several parts.

### 5.1. Lid assembly

### 5.2. Instructions for assembling the cell for the **side-by-side** arrangement of electrodes.


### 5.3. Instructions for assembling the cell for the **face-to-face** arrangement of electrodes.

These instructions can also be viewed as videos on our website [el-cell.com](http://el-cell.com).

## 5.1 Lid assembly:




1. Turn the lid upside down.  
Insert the foil sealing into the lid.



2. Insert the window.



3. Insert the window thrust ring.



4. Attach the window flange and screw in the screws using the torque screwdriver 0.2 Nm.




## 5.2 Assembling the cell with a side-by-side sample holder:




1. Assemble the lid as described above. Then insert the lower piston into the inner sleeve.




2. Use the loading tweezers to push the piston into the sleeve.



3. Place the assembly into the cell base.



4. Insert the sample holder into the loading tweezers. Squeeze the tweezers together to lift the sample holder bracket. Then attach the first electrode strip to the sample holder.



5. Turn the sample holder around and place the second electrode strip.



6. Put the R electrode disc, lithium metal in this case, on top of the piston.



7. Put the separator on top. The separator must not be more than 10 mm in diameter. A smaller separator may be useful to prevent the sample holder from coming into contact with the electrolyte.



8. Use the loading tweezers again to press the cell stack further into the piston.



9. Attach the sample holder to the assembly. The markings 1 and 2 on the sample holder and cell base represent the electrodes. Use them as an aid for correct placement.



10. Insert the second foil seal.



11. Drop the electrolyte onto the separator/electrodes. Use only the amount of electrolyte necessary to impregnate the porous material. Avoid excess electrolytes. Typically, less than 50  $\mu$ l of electrolyte is suitable.

Note that electrolyte will be squeezed out of the separator when the cell lid is attached.



12. Use the straight marking on the lid to align it to the flat side of the cell base and attach it.



13. Finally fix the lid by using the Torx screws and the hexagon screwdriver 2.5mm.



14. The PAT-Cell-Opto-10 is now ready for testing.

### 5.3 Assembling the cell with a face-to-face sample holder:



1. Follow the lid assembly steps. Then insert the lower piston into the inner sleeve.




2. Use the loading tweezers to push the piston into the sleeve.



3. Place the assembly into the cell base.




4. Put the R electrode disc, lithium metal in this case, on top of the piston.



5. Put the separator on top. The separator must not be more than 10 mm in diameter. A smaller separator may be useful to prevent the sample holder from coming into contact with the electrolyte.



6. Use the loading tweezers again to push the cell stack further down.



7. Drop the electrolyte onto the separator. Use only the amount of electrolyte necessary to impregnate the porous material. Avoid excess electrolytes. Typically, less than 50  $\mu$ l of electrolyte is suitable.

Note that electrolyte will be squeezed out of the separator when the cell lid is attached.



8. Put the self-standing electrode on top.



9. Push down the electrode stack using the loading tweezers.




10. Insert the contact ring.




11. Put the holed current collector on top. Choose the hole size depending on the magnification used by the microscope.




12. Insert the second foil seal.



13. Put on the cell lid. Make sure that the flat side of the cell base and the straight line on the lid are parallel.



14. Finally, screw the cover on with the 2.5 mm hexagon screwdriver.



15. The PAT-Cell-Opto-10 is now ready for testing.

## 6 Connecting the test cell

The PAT-Cell-Opto-10 is connected via the cableless PAT socket. It can be directly plugged into any EL-Cell PAT-Tester or PAT docking station.



## 7 Disassembly and Cleaning

Right after use, disassemble the cell in the reverse order of assembly. All chemicals used have to be disposed of properly. All wetted parts are to be cleaned with water and/or other appropriate solvents. All parts are to be dried immediately after cleaning at 80°C. All cell parts made of PPS are to be dried immediately after cleaning overnight at 120°C.

**NOTE:** Leaving cell parts in contact with the ambient atmosphere while still being wetted with electrolyte may result in severe corrosion.

## 8 Unpacking

Check the contents of the packages against the list given below to verify that you have received all of the components. Contact the factory if anything is missing or damaged.

**NOTE:** Damaged shipments must remain with the original packaging for freight company inspection.



### List of Components:

1. **PAT-Cell-Opto-10** test cell with flat socket, fully equipped for use in both 2-electrode and 3-electrode (reference) configuration
2. **Windows seal (10 pcs.)** [ECC1-05-0016-B/X](#)
3. **Separator 10.0 mm x 0.26 mm, GF/A (10 pcs.)** [ECC1-01-0012-R/X](#)
4. **PE Sealing foil (10 pcs.)** [ECC1-05-0032-A/X](#)
5. **Loading tweezer** [ECC1-09-2010-A](#)
6. **Tweezer antiaacid / stainless (pointed)** [WZG9020](#)
7. **Torque screwdriver, 0.2 Nm, cross handle** [WZG9021](#)
8. **Inbus bit 1.5 mm x ¼ inch** [WZG9046](#)
9. **Hexagon screwdriver 2.5 mm, cross handle** [WZG9047](#)
10. **Sample holder (side-by-side)** [ECC1-05-0039-A](#)
11. **Contact ring (face-to-face)** [ECC1-05-0039-C](#)
12. **Contact disc (50 fold)** [ECC1-05-0042-C](#)
13. **Contact disc (100 fold)** [ECC1-05-0042-A](#)



14. **Contact disc (200 fold)** [ECC1-05-0042-B](#)

## 9 Components

### PAT-Cell-Opto-10



## Lid (OPTO-10) 18 mm, assy



## 10 Technical Support

Technical support for this product is exclusively provided by EL-Cell GmbH.

### **EL-Cell GmbH**

Tempowerkring 8

21079 Hamburg - Germany

**phone:** +49 40 79012-737

**fax:** +49 40 79012-736

**e-mail:** [support@el-cell.com](mailto:support@el-cell.com)

**web:** [www.el-cell.com](http://www.el-cell.com)

## 11 Warranty

For a period of one year from the date of shipment, EL-Cell GmbH (hereinafter Seller) warrants the goods to be free from defects in material and workmanship to the original purchaser. During the warranty period, Seller agrees to repair or replace defective and/or nonconforming goods or parts without charge for material or labor, or, at the Seller's option, demand return of the goods and tender repayment of the price. The buyer's exclusive remedy is repair or replacement of defective and nonconforming goods, or, at the Seller's option, the repayment of the price.

Seller excludes and disclaims any liability for lost profits, personal injury, interruption of service, or consequential incidental or special damages arising out of, resulting from, or relating in any manner to these goods.

This Limited Warranty does not cover defects, damage, or nonconformity resulting from abuse, misuse, neglect, lack of reasonable care, modification, or the attachment of improper devices to the goods. This Limited Warranty does not cover expendable items. This warranty is void when repairs are performed by a non-authorized person or service center. At the Seller's option, repairs or replacements will be made on-site or at the factory. If repairs or replacements are to be made at the factory, the Buyer shall return the goods prepaid and bear all the risks of loss until delivered to the factory. If Seller returns the goods, they will be delivered prepaid and Seller will bear all risks of loss until delivery to Buyer. Buyer and Seller agree that this Limited Warranty shall be governed by and construed by the laws of Germany.

The warranties contained in this agreement are in lieu of all other warranties expressed or implied, including the warranties of merchantability and fitness for a particular purpose.

This Limited Warranty supersedes all prior proposals or representations oral or written and constitutes the entire understanding regarding the warranties made by Seller to Buyer. This Limited Warranty may not be expanded or modified except in writing signed by the parties hereto.