

# User Manual

Release 1.12

## PAT-Terminal-1



## Disclaimer

EL-Cell GmbH makes no assurances or warranties with respect to this manual and, to the extent permitted by law, limits its liability for violation of any implied warranty to the substitution of this manual for another. In addition, EL-Cell GmbH reserves the right to revise this publication at any time without notice to anyone.

The information provided in this documentation includes general descriptions and/or technical characteristics regarding the performance of the equipment described herein. This documentation cannot serve as a proper evaluation of the suitability or reliability of the equipment for any specific application by any user and should not be relied upon as a substitute for such evaluation. It is the responsibility of each such user or installer to conduct an appropriate and complete risk assessment, evaluation and testing of the equipment with respect to their specific application. EL-Cell GmbH cannot be held responsible or liable for any misuse of the information contained herein.

All relevant state, regional and local safety regulations must always be complied with when installing and using this device. For safety reasons and to ensure compliance with the documented system data, only the manufacturer is authorized to perform repairs on components.

Disregarding this information may result in injury or damage to the equipment.

All rights reserved. No part of this publication may be reproduced, processed, or transmitted in any form, including photocopying, recording, or any other electronic or mechanical process, without the written permission of the publisher. Requests for permission must be made in writing to the publisher at the address written below.

## Manufacturer and customer service

EL-Cell GmbH  
Tempowerkring 8  
21079 Hamburg – Germany

**Telephone:** +49 40 79012-734

**Telefax:** +49 40 79012-736

**Email:** [info@el-cell.com](mailto:info@el-cell.com)

**Website:** [el-cell.com](http://el-cell.com)

## Technical support

**Telephone:** +49 40 79012-734

**Email:** [support@el-cell.com](mailto:support@el-cell.com)

**Website:** [el-cell.com/support/technical-support/](http://el-cell.com/support/technical-support/)

Please always quote the serial number on the nameplate when making customer service inquiries.

## Shipping address for repairs

EL-Cell GmbH  
Tempowerkring 8  
21079 Hamburg - Germany

Please be sure to contact our customer service department before making a return. We will not open or process shipments without a completed decontamination report or RMA.

## Table of contents

<b>1</b>	Preamble.....	6
1.1	Purpose and target audience.....	6
1.2	Storage instructions.....	6
1.3	Obtaining documents and information.....	7
<b>2</b>	For your safety.....	8
2.1	Explanation of the safety instructions.....	8
2.1.1	Terms used.....	8
2.1.2	Symbols used.....	8
2.2	Product safety and hazards.....	9
2.2.1	Requirements for operating personnel.....	9
2.3	Owner's responsibility.....	9
2.4	Intended use.....	10
2.5	Modifications and conversions.....	10
2.6	Behavior in the event of malfunctions and irregularities.....	10
2.7	Switch off device in emergency.....	10
<b>3</b>	Description and design.....	11
3.1	Description.....	11
3.2	Design and connections.....	11
3.2.1	Connections.....	12
3.2.2	Pin assignment D-Sub connector.....	13
3.3	Technical data.....	14
3.3.1	Dimensions.....	14
3.3.2	General device specifications.....	14
3.3.3	Test channel performance data.....	15
3.4	Applied guidelines and standards.....	17
	The product described is in conformity with the following harmonized standards:.....	17
3.5	Checking for completeness and transport damage.....	20
3.6	Scope of delivery.....	20
3.7	Recycling of the packaging material.....	20

3.8	Storage after delivery .....	20
4	Start-up .....	21
4.1	Operation as a stand-alone device .....	21
4.1.1	Operation outside the glove box .....	21
4.1.2	Operation inside the glove box .....	21
4.1.3	Performing an impedance test.....	22
4.2	Operation as test channel in a PAT-Tester-x setup .....	22
4.3	Switching on and off.....	22
4.4	Inserting and removing test cells .....	23
4.4.1	Inserting a PAT series test cell .....	23
4.4.2	Removing a PAT series test cell .....	23
4.4.3	Using other cell types .....	23
5	LC-Display.....	24
6	Malfunctions, warning and error messages .....	25
6.1	Explanation of the visual signals.....	25
6.1.1	LED Front panel .....	25
6.1.2	LED PAT socket .....	26
6.1.3	Troubleshooting .....	26
7	Maintenance and repair .....	27
7.1	Cleaning.....	27
8	Storage and disposal .....	28
8.1	Storage.....	28
8.2	Disposal .....	28
9	Warranty .....	29

# 1 Preamble

## 1.1 Purpose and target audience

This manual describes the structure, function, operation and maintenance of the PAT-Terminal-1. It is intended for the end users of the device. An end user can be described as any person who interacts directly with the PAT-Terminal-1. The term "end user" usually includes laboratory personnel who have been specifically trained to operate this instrument and are familiar with all the precautions required to work in the laboratory.

Only an authorized and properly qualified and experienced person 18 years of age or older may use the PAT-Terminal-1 machine, who:

- has read and understood these installation and operating instructions
- is familiar with the installation and operation of this or a similar device
- is aware of all possible dangers and acts accordingly

## 1.2 Storage instructions

Make sure you have read and understood the complete instructions and all safety information before using this product. Failure to follow these instructions may result in minor or serious injury.

Follow all instructions. This will prevent accidents that could result in property damage or injury. Keep all safety information and instructions for future reference and pass them on to subsequent users of the product.

The manufacturer is not liable for property damage or injuries resulting from incorrect handling or failure to comply with the safety instructions. In such cases, the warranty becomes void.

## 1.3 Obtaining documents and information

A current version of the documentation is available on the following website:

<https://el-cell.com/support/manuals/>

Alternatively, you can scan this QR code,  
to access the website:



## 2 For your safety



### 2.1 Explanation of the safety instructions

In this manual, certain recurring terms and symbols are used to warn you of hazards or to give you instructions that are important to prevent injury and damage. It is essential that you read and follow these notes in order to avoid accidents and damage. These terms and symbols are explained below.

#### 2.1.1 Terms used

NOTICE
„Notice" indicates an important or useful additional information

#### 2.1.2 Symbols used

Warning symbols (warn of a hazard)	Mandatory sign (prescribes an action)
 Risk of electric shock	 Pull out power plug



## 2.2 Product safety and hazards

The device described is technically mature, is manufactured using high-quality materials and is tested at the factory before delivery. It complies with the state of the art and the recognized safety regulations.

### 2.2.1 Requirements for operating personnel

The device may only be operated and maintained by persons of legal minimum age who have been instructed in its use.

## 2.3 Owner's responsibility

The owner of the device

- is responsible for the proper condition of the device and for ensuring that it is operated in accordance with its intended use;
- is responsible for ensuring that persons who are to operate or maintain the device are technically qualified to do so, are instructed on the device and are familiarized with these operating instructions;
- must be familiar with the applicable rules, regulations and occupational health and safety provisions and train the personnel accordingly;
- is responsible for ensuring that unauthorized persons do not have access to the device;
- is responsible for ensuring that personal protective equipment is worn by the operating personnel, e.g. work clothing, safety shoes, protective gloves.

## 2.4 Intended use

The PAT-Terminal-1 is designed for charging and discharging, as well as for performing functional tests on electrochemical battery test cells. The device can be operated inside a glove box under a protective atmosphere. It can be used independently or in conjunction with a PAT-Tester-x. Other uses can lead to hazards and damage.

### NOTICE

The PAT-Terminal-1 may only be used with test cells of the PAT series. Other cell types must be connected using the appropriate adapters available from the manufacturer.

## 2.5 Modifications and conversions

The PAT-Terminal-1 must not be modified or altered without authorization. No parts may be added or installed that are not approved by the manufacturer.

Unauthorized conversions or modifications will result in the device no longer being allowed to be operated.

The manufacturer is not liable for any damage, danger or injury resulting from unauthorized conversions or modifications or from non-compliance with the regulations in this manual.

## 2.6 Behavior in the event of malfunctions and irregularities

The PAT-Terminal-1 may only be operated in perfect condition. If you, as the operator, notice irregularities, malfunctions or damage, take the device out of operation immediately and inform your supervisor.

## 2.7 Switch off device in emergency

In case of emergency, unplug the power supply unit to disconnect the device from the mains supply.

## 3 Description and design


### 3.1 Description

The PAT-Terminal-1 is a single-channel, fully equipped potentiostat / galvanostat and impedance analyzer. It has an LC display for displaying various parameters of the test cell used. The PAT-Terminal-1 can be operated as a stand-alone unit in a glove box to perform functional tests on test cells and to facilitate cell sensor adjustment. Once connected to a PAT Controller-8, it can be used as a normal test channel for full electrochemical measurements.

### 3.2 Design and connections

Front side


Back side



1	Eject button	4	Status LED, device
2	Function key	5	Status LED, PAT socket
3	LC-Display		


### 3.2.1 Connections

Front side




1. **PAT-Channel connector (USB 2 Type C):** For connecting another PAT-Channel-1
2. **AUX ports (USB 2.0):** For future applications

Back side




3. **PAT-Controller connector (USB 2 Type C):** Used to supply power via the supplied USB power adapter or to connect to a PAT-Controller.

Top side




4. **D-Sub connector:** This connector is used to connect a battery test cell or docking station via D-Sub cable.
5. **USB 2.0 port:** Used for power supply and data transfer of external devices.
6. **PAT socket:** This socket is used to connect a PAT series battery test cell or PAT-Adapter.

### 3.2.2 Pin assignment D-Sub connector



### 3.3 Technical data



#### 3.3.1 Dimensions

Weight	1.3 kg (without inserted test cell)
Height	97 mm
Length	164 mm
Width	141 mm

#### 3.3.2 General device specifications

Device name	PAT-Terminal-1
Type	Potentiostat/Galvanostat/Impedance analyzer
Operating temperature range	-20 to +40 °C
Connections	<ul style="list-style-type: none"> <li>• D-Sub connector for actively shielded cell cables, I2C bus and analog signals</li> <li>• USB 2.0 connectors</li> <li>• USB 2 Type C connector for power supply and connection to a PAT controller-8</li> </ul>

### 3.3.3 Test channel performance data

General	
Control Voltage	-7 to +7 V
Compliance Voltage	-8 V to +8 V (no load)
Current	±100 mA
Cell and electrode connections	3 electrodes plus sense wires, Connection matrix
ADC	2x24 Bit
DAC	1x18 Bit
Bandwidth ranges	500 kHz 50 kHz 5 kHz
Slew rate	2.5 V/μs
Sampling interval (rate)	1 ms
Input Impedance	>100 MΩ    20 pF
Computer Interface	1 GBit Ethernet Runs standalone Multiuser
Voltage	
Acquisition voltages	Full cell voltage Both half cell voltages Auxiliary voltage
Measurement Accuracy	±0.02% FSR (Full Scale Range)
Measurement Noise floor	30 μV peak-peak typical
Control Resolution	57 μV (18 Bit) EIS amplitude: 3 μV (additional 16 Bit DAC for EIS)
Current	
Current Ranges	±100 mA ±10 mA ±1 mA ±100 μA Auto Range
Measurement Noise floor	<1 μA @ 100mA <100 nA @ 10mA

	<10 nA @ 1mA <1 nA @ 100µA
Measurement Accuracy	±0.05% FSR
Control Resolution	1 nA min. (18 Bit)
<b>Impedance (each channel)</b>	
Frequency range	100 µHz to 100 kHz
Impedance mode	PEIS and GEIS (simultaneous measurement of full- and half-cell impedances)
Impedance range	1 mΩ to 100 MΩ
EIS quality indicator	SFDR (Spurious Free Dynamic Range)
EIS drift correction	yes
EIS adaptive amplitude	yes
<b>Other</b>	
Additional Measurement (each channel)	Digital (I <sup>2</sup> C) sensor signal, e.g. for cell temperature Analog sensor signal, e.g. for gas pressure
Calibration	Fully automatic self-calibration with internal voltage reference and three internal calibration cells
Cell Identification	PAT-Button with unique serial number stored in EEPROM



### 3.4 Applied guidelines and standards

The product described is in conformity with the following harmonized standards:

<b>EN 61010-1:2010</b>	Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte – Teil 1: Allgemeine Anforderungen (DIN EN 61010-1, VDE 0411-1:2011-07)
	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements (IEC 61010-1:2010 + Cor. :2011)
<b>EN 61010-2-201:2014</b>	Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte - Teil 2-201: Besondere Anforderungen für Steuer- und Regelgeräte (DIN EN 61010-2-201:2014, VDE 0411-2-201:2014-01)
	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-201: Particular requirements for control equipment (IEC 61010-2-201:2013)
<b>EN 61010-2-010:2015-05</b>	Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte - Teil 2-010: Besondere Anforderungen an Laborgeräte für das Erhitzen von Stoffen (DIN EN 61010-2-010:2014; VDE 0411-2-010:2015-05)
	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-201: Particular requirements for control equipment (IEC 61010-2-201:2013)
<b>EN 61326-1:2013</b>	Elektrische Mess-, Steuer-, Regel- und Laborgeräte - EMV-Anforderungen - Teil 1: Allgemeine Anforderungen (DIN EN 61326-1:2013-07, VDE 0843-20-1:2013-07)
	EMC requirements - Part 2-3: Particular requirements - Test configuration, operational conditions and performance criteria for transducers with integrated or remote signal conditioning

	(IEC 61326-2-3:2012)
<b>EN 61326-2-3:2013-07</b>	Elektrische Mess-, Steuer-, Regel- und Laborgeräte - EMV-Anforderungen - Teil 2-3: Besondere Anforderungen - Prüfanordnung, Betriebsbedingungen und Leistungsmerkmale für Messgrößenumformer mit integrierter oder abgesetzter Signalaufbereitung  (DIN EN 61326-2-3:2013-07, VDE 0843-20-2-3:2013-07)
	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-3: Particular requirements - Test configuration, operational conditions and performance criteria for transducers with integrated or remote signal conditioning  (IEC 61326-2-3:2012)
<b>EN 50581: 2013-02</b>	Technische Dokumentation zur Beurteilung von Elektro- und Elektronikgeräten hinsichtlich der Beschränkung gefährlicher Stoffe  (DIN EN 50581; VDE 0042-12:2013-02)
	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances



electrochemical test equipment

## EU Declaration of Conformity

Manufacturer's name and address: EL-Cell GmbH  
Tempowerkring 8  
21079 Hamburg  
Germany

Product: PAT-Terminal-1

**The designated product is in conformity with the**

- Low Voltage Directive (LDV) 2014/35/EU
- Electromagnetic Compatibility Directive (CEM) 2014/30/EU
- Restriction of Hazardous Substance Directive (RoHS) 2011/65/EU

**and the following harmonised standards:**

- Safety: IEC 61010-1
- EMC: IEC 61326

Emissions

EN 55011: Conducted Class B  
EN 55011: Radiated Class A  
EN 61000-3-2: Harmonic Current

Immunity

IEC 61000-4-3: EM field  
IEC 61000-4-4: Burst  
IEC 61000-4-5: Surge  
IEC 61000-4-6: Conducted RF  
IEC 61000-4-8: Magnetic Field  
IEC 61000-4-11: Voltage Dip/Short Interruptions

Hamburg, 25.05.2023

Jan Kahrs, CEO

This declaration certifies compliance with the above mentioned directives but does not include a property assurance.  
The safety note given in the product documentation which are part of the supply, must be observed.

## 3.5 Checking for completeness and transport damage

- Check the completeness of the scope of delivery against the delivery note.
- Check the device for damage.

If you find deviations from the scope of delivery or damage, please notify the carrier and the manufacturer.

## 3.6 Scope of delivery

- PAT-Terminal-1
- USB 2.0 Typ C Kabel, 1m, [ELT9798](#)
- Power cable, 0.75m, [ELT1053](#)
- USB-C Power supply 65W, [ELT1045](#)

## 3.7 Recycling of the packaging material

Dispose of the packaging material (cardboard, foil) in accordance with the legal regulations applicable in your country.

## 3.8 Storage after delivery

If the device is to be stored initially after delivery, please observe the storage conditions.

## 4 Start-up

### 4.1 Operation as a stand-alone device

#### 4.1.1 Operation outside the glove box

1. Place the PAT-Terminal-1 on a horizontal, dry and clean surface.
2. Connect the supplied power supply unit to the PAT-Terminal-1 and then connect the power supply unit to the mains.
3. You can now insert a battery test cell into the PAT socket. The electrical values of the test cell appear on the display.

#### NOTICE

After removing the battery test cell, values continue to be shown on the display until a new cell is connected to the device. This is an intended behavior.

#### 4.1.2 Operation inside the glove box

1. Insert the PAT-Terminal-1 into the glove box (transport through vacuum permitted).
2. Now connect the device to the supplied power supply outside the glove box. To do this, the supplied USB feedthrough-4 must be installed on your glove box.
3. First connect the PAT-Terminal-1 to the cable feedthrough using the USB cable supplied.
4. Then connect the USB power supply to the cable feedthrough and to the mains.
5. You can now insert a battery test cell into the PAT socket. The electrical values of the test cell appear on the display.

#### NOTICE

If the PAT-Terminal-1 is not supplied with power even though the cables are connected, disconnect the USB cable in the glove box from the feedthrough, turn it 180 degrees and plug it back in.

### 4.1.3 Performing an impedance test

By pressing the function key next to the display, an impedance measurement (PEIS, 1 KHz, 10 mV) is triggered and the cell impedances Z12, Z1 and Z2 are output on the display.

## 4.2 Operation as test channel in a PAT-Tester-x setup

### Precondition

To operate the PAT-Terminal-1 as a test channel in a PAT-Tester-x setup, you need a PAT-Controller connected to an Ethernet Local Area Network (LAN). EL-Software Server (the server application of the EL-Software) must be installed on the LAN server and EL-Software Client (the client component of the EL-Software) must be installed on at least one client PC. The installation instructions for the PAT-Controller and the EL software can be found in separate manuals.

### IMPORTANT NOTICE

Please make sure that the firmware of the PAT-Controller is up to date before connecting the PAT-Terminal to it. For more information on how to perform a firmware upgrade, please refer to the EL-Software manual.

1. Place the PAT-Terminal-1 on a flat, dry and clean surface.
2. Connect the PAT-Terminal-1 to the PAT-Controller-8 using the supplied USB cable. The cable connects the USB port on the PAT-Terminal-1 (labeled "PAT-Controller") to one of the 8 USB ports (labeled "Channel") on the PAT-Controller.
3. Plug a PAT test cell into the PAT socket. Alternatively, you can also connect another cell type or a PAT docking station via the D-Sub connector of the PAT-Channel and a special cable. This option can be used, for example, to operate a test cell at extreme temperatures in a temperature chamber. An adapter is optionally available for testing button cells in the PAT socket.


### 4.3 Switching on and off

The PAT-Terminal-1 switches on as soon as it is connected to the power supply. To switch off, it is sufficient to disconnect the device from the power supply.

## 4.4 Inserting and removing test cells


### 4.4.1 Inserting a PAT series test cell

Insert the cell into the PAT socket until the cell audibly clicks into place. This ensures that the electrical contact is properly established.



### 4.4.2 Removing a PAT series test cell

To remove a cell, press the eject button next to the PAT socket to release the lock (1) and then remove the cell (2).



### 4.4.3 Using other cell types

Other cell types can be connected to the PAT socket via an adapter or to the D-Sub input.



## 5 LC-Display



Display item	Meaning
Cell ID (in this case F0001)	Cell ID of the test cell used
Force	The force acting on the cell stack in Newtons (if sensor is present).
Press	The gas pressure in the test cell in mbar (if sensor is present)
Temp	The temperature in the test cell in °C (if sensor is present)
V12, V1R, V2R	The voltages of the full cell and the two half cells in Volts
Z12, Z1, Z2	The impedances of the full cell and the two half cells at 1 kHz in Ohms (only after pressing the function key)



## 6 Malfunctions, warning and error messages

WARNING	
	After removing covers, electrically live parts may be exposed. You may suffer an electric shock if you touch these parts.
	Disconnect the power plug before removing any covers. Only qualified electricians may work on the electrical equipment of the devices.

### 6.1 Explanation of the visual signals

#### 6.1.1 LED Front panel

LED	Signal	Meaning
Front panel	Red, lighting interval 3 sec	Initialising Mainboard
Front panel	Red, permanent	Mainboard Error
Front panel	Green, permanent	Mainboard operational
Front panel	Blue, permanent	PAT-Controller is waiting for assignment of an IP address from the DHCP server <b>or</b> Establishing USB connection <b>or</b> Performing firmware update (Mainboard+Channelboard, can take up to 5 min)

## 6.1.2 LED PAT socket

Signal	Bedeutung
Red, flashes every three seconds	Channelboard is initialized
Red, permanent	Channelboard error
Green, permanent	Channelboard works normally

## 6.1.3 Troubleshooting

### Problem:

The PAT-Terminal-1 displays an empty display after connecting to a PAT-Controller-8 and does not react anymore.

### Possible Solution:

The firmware version on the PAT-Controller-8 is too old and must be updated. This is the case if the installed EL-Software release is older than version 1.1.43.

Proceed as follows:

1. Update EL-Software to the latest version. Detailed information about the update process can be found in the EL-Software Manual.
2. Then update the firmware of the PAT-Controller-8.
3. Reconnect the PAT-Terminal-1 to the PAT-Controller-8. It will now receive the latest firmware and should work properly afterwards.

## 7 Maintenance and repair

### 7.1 Cleaning

Wipe the PAT Terminal-1 with a moist cloth. Do not use aggressive chemicals for cleaning. The metal surfaces can be cleaned with commercially available stainless steel cleaning agents. If rust spots occur due to contamination, the affected areas must be cleaned and polished immediately.

Protect the PAT-Terminal-1 from dust and splash water.

## 8 Storage and disposal

### 8.1 Storage

The device may only be stored under the following conditions:

- dry and in a closed, dust-free room
- frost-free
- disconnected from the power supply

### 8.2 Disposal

The device must not be disposed of in normal household waste. Observe the applicable legal regulations. Please contact your dealer or the manufacturer for disposal.

## 9 Warranty

For a period of one year from the date of shipment, EL-Cell GmbH (hereinafter Seller) warrants the goods to be free from defect in material and workmanship to the original purchaser. During the warranty period, Seller agrees to repair or replace defective and/or nonconforming goods or parts without charge for material or labor, or, at the Seller's option, demand return of the goods and tender repayment of the price. Buyer's exclusive remedy is repair or replacement of defective and nonconforming goods, or, at Seller's option, the repayment of the price.

Seller excludes and disclaims any liability for lost profits, personal injury, interruption of service, or for consequential incidental or special damages arising out of, resulting from, or relating in any manner to these goods.

This Limited Warranty does not cover defects, damage, or nonconformity resulting from abuse, misuse, neglect, lack of reasonable care, modification, or the attachment of improper devices to the goods. This Limited Warranty does not cover expendable items. This warranty is void when repairs are performed by a non-authorized person or service center. At Seller's option, repairs or replacements will be made on site or at the factory. If repairs or replacements are to be made at the factory, Buyer shall return the goods prepaid and bear all the risks of loss until delivered to the factory. If Seller returns the goods, they will be delivered prepaid and Seller will bear all risks of loss until delivery to Buyer. Buyer and Seller agree that this Limited Warranty shall be governed by and construed in accordance with the laws of Germany.

The warranties contained in this agreement are in lieu of all other warranties expressed or implied, including the warranties of merchantability and fitness for a particular purpose.

This Limited Warranty supersedes all prior proposals or representations oral or written and constitutes the entire understanding regarding the warranties made by Seller to Buyer. This Limited Warranty may not be expanded or modified except in writing signed by the parties hereto.